

Capsule Chart Xpress™

PRODUCT BRIEF



Enhances clinical decision-making

- Faster entry of vital signs into the medical record
- Accurate charting of vitals and patient-specific information
- Identify patients potentially at risk using the optional Capsule Early Warning Scoring System

Increases time for direct care

- Simplifies established clinical work flows
- Vital signs are filed to the patient record at the point-of-care

Cost-effective

- Adds connectivity and documentation to all leading spot check monitors

Reduce Charting Time from Hours to Seconds

Are you capturing and manually documenting vital signs in med-surg? Do they often end up on exam gloves, alcohol wipes and scraps of paper? And does it take hours to transcribe or key those vitals into the electronic medical record?

With Capsule Chart Xpress™, periodic vital signs enter the patient's record electronically, before the clinician leaves the bedside. Utilizing the Capsule Neuron™ 2 mobile clinical computer, Chart Xpress connects to your existing spot check monitors and wirelessly sends validated vital signs to the electronic medical record. Vitals are reviewed and validated at the point-of-care, so clinicians don't have to log in to another workstation or carry another handheld device.

Designed with nurses for nurses, Chart Xpress is intuitive, flexible and easy to use. It can be configured to any clinical workflow, and can accept additional vital signs parameters unique to the care environment.

The result? Chart Xpress eliminates time-consuming, error-prone manual entry of vitals during or after rounds, literally reducing charting time from hours to seconds. The care team can devote more time to direct care activities and spend less time in charting.



Capsule Chart Xpress is a product of Qualcomm Life, Inc.

Capsule ChartXpress Workflow

Vital Signs Documentation in 6 Simple Steps

- 1 Login**
- 2 Identify Patient**

405	ID : GT45612
JOHNSON, Kenneth	
Gender	Male
Birth	Jul 5, 1978
Age	34
ID	GT45612

I have confirmed this patient's identity using two acceptable forms of patient identifiers
- 3 Capture Vitals**

Temp	37.2 °C	Oral	Tympanic	Axillary	Rectal
Pulse	78 bpm	Monitor	Apical	Right Radial	Left Radial
Resp	24	Observed			
BP	170/80 mm Hg	Lying	Sitting	Standing	R Arm
MAP	110 mm Hg				L Arm
O2Sat	98 %	Room air	Nasal cannula	Simple mask	Venturi
- 4 Document**

AVPU Scale	Alert	Voice	Pain	Unrespons...
Fall Precaution	SRx4	SRx2	Bed Low	Wheels loc...
Weight	154	lbs		

Comments: Text entered here
- 5 Review**

Temp	37.2 °C	Oral	AVPU	Pain
Pulse	78 bpm	Right Radial	Fall precaution	Bed low
Resp	24 rr min	Observed	Weight	154 lbs
BP	170/80 mm Hg	Standing		
MAP	110 mm Hg	L Arm		
O2 Sat	98 %	Venturi		
- 6 Send**

Vitals Submitted on 05 Feb 2013 17:41:07

CONTINUE SAME PATIENT | SELECT NEW PATIENT

For more information, contact us:

In North America

+978.482.2300
 qcl.capsule.support@qualcommmlife.com

International Offices

+33 1 84 17 12 68
 qcl.international@qualcommmlife.com

Visit qualcommmlife.com

Capsule Chart Xpress™ is an application included in Capsule's Medical Device Information System (MDIS)—a platform of products designed to increase efficiency and improve patient care delivery through the systematic, timely capture of patient vital signs at the point-of-care, and the delivery of captured data to hospital electronic medical records, other information systems and third party mobile solutions.